# **UNCLASSIFIED**

# AD NUMBER AD399956 CLASSIFICATION CHANGES TO: unclassified FROM: confidential LIMITATION CHANGES

# TO:

Approved for public release; distribution is unlimited.

# FROM:

Distribution authorized to DoD only; Foreign Government Information; FEB 1970. Other requests shall be referred to British Embassy, 3100 Massachusetts Avenue, NW, Washington, DC 20008.

# AUTHORITY

DSTL ltr dtd 12 Dec 2006; DSTL ltr dtd 12 Dec 2006

ERDE TN 10

CONFIDENTIAL/DISCREET (but cleared for Canada)



MINISTRY OF TECHNOLOGY

Decl OADR

This material centains information affecting the national defense of the United States within the mosning of the Espionage Laws (Title 18, U.S.C., sections 783 and 794), the transmission or reveletion of which in any manner to an unauthorized person is prohibited by law.

# XPLOSIVES RESEARCH **DEVELOPMENT ESTABLISHMENT**

TECHNICAL NOTE No. 10

PICATINNY ARSENAL

EXCLUDED FROM AUTOMOTIO REGRADIES; DOD PIE 5200.10 DOES NOT APPLE

The Spectrophotometric Determination of Resorcinol in Colloidal Propellants (C)

# 20090108 005

# F.I.H. Tunstall

### THIS DOCUMENT IS THE PROPERTY OF HER BRITANNIC MAJESTY'S GOVERNMENT

It is issued for the information of such persons only as need to know its contents in the course of their official duties. Any person other than the authorised holder, upon obtaining possession of this document by finding or otherwise, should forward it together with his name and address in a closed envelope to:

SECURITY BRANCH, MINISTRY OF TECHNOLOGY. ST. GILES COURT, ST. GILES HIGH STREET, LONDON, W.C.2.

Letter postage need not be prepaid, other postage will be refunded. ALL PERSONS ARE HEREBY WARNED THAT THE UNAUTHORISED RETENTION OR DESTRUCTION OF THIS DOCUMENT IS AN OFFENCE AGAINST THE OFFICIAL SECRETS ACT.

FOR OVERSEAS RELEASE CONDITIONS SEE INSIDE COVER

WALTHAM ABBEY **ESSEX** 

CONFIDENTIAL/DISCREET

(but cleared for Canada)

+# 115572

956668-07

# MINISTRY OF TECHNOLOGY

EXPLOSIVES RESEARCH AND DEVELOPMENT ESTABLISHMENT, G. Brit.

Technical Note No. 10

8.1

6 | February 1970)

The Spectrophotometric Determination of Resorcinol in Colloidal Propellants (C)

by

F.I.H. Tunstally

### SUMMARY

The method for determination of resorcinol which is based on spectrophotometric measurement of copper nitrosoresorcinate has been modified to improve its sensitivity and reliability.

Resorcinols. Tuntall F. 2 Double base rocket propellants. 16x. Brit- Resorcinols. 2. Spectrophotometry. 16x. Brit- Double base rocket propellants.

CONFIDENTIAL/DISCREET (but cleared for Canada)

V. To the said

Further copies of this technical note can be obtained from Mintech Reports Centre, Station Square House, St. Mary Cray, Orpington, Kent BR5 3RE

CONFIDENTIAL/DISCREET
(but cleared for Canada)

# CONTENTS

|    |                | Page No |
|----|----------------|---------|
| 1. | Introduction   | 1       |
| 2. | Experimental   | 2       |
| 3. | Procedure      | 3       |
| 4. | Results        | 4       |
| 5. | Conclusions    | 4       |
| 6. | References     | 5       |
|    | Figures 1 to 4 |         |

Reference: WAC/205/08

# 1. INTRODUCTION

To determine resorcinol by the method given in Reference 1, the propellant is first extracted with dichloromethane. Resorcinol is then separated from other soluble ingredients including NG, triacetin, and 2-NDPA by washing with water. Treatment of the aqueous extract with sodium nitrite and dilute acetic acid produces 2:4-dinitrosoresorcinol (Eqn. 1) and addition of copper sulphate solution produces a yellow copper dinitrosoresorcinate complex (Eqn. 2). This exhibits two absorption maxima at 290 and 340 mm and measurement at 340 mm enables resorcinol to be determined. The structure of the complex is indefinite and the reaction probably proceeds through several equilibria as follows:

$$\begin{array}{c}
OH \\
NaNO_2 \\
HOAc
\end{array}$$

$$\begin{array}{c}
NaNO_2 \\
HOAc
\end{array}$$

$$\begin{array}{c}
OH \\
NO
\end{array}$$

$$\begin{array}{c}
OH \\
NO$$

$$\begin{array}{c}
OH \\
NO
\end{array}$$

$$\begin{array}{c}
OH \\
NO
\end{array}$$

$$\begin{array}{c}
OH \\
NO$$

$$\begin{array}{c}
OH \\
NO
\end{array}$$

$$\begin{array}{c}
OH \\
NO$$

$$\begin{array}$$

The considerable water solubility of resorcinol enables it to be separated from almost all other propellant ingredients which might interfere with the estimation and the method can be used satisfactorily to analyse degraded propellants from ageing trials. Since the work was completed several authors have described chromatographic procedures for the determination of resorcinol in CMCDB propellant compositions.<sup>2,3,4</sup> Propellant samples from compatibility tests which are mixed with rubber and other materials have also been analysed successfully, although in such work it is necessary to prove in each case that the method is unaffected by this contamination. It is sometimes preferable to use diethyl ether rather than dichloromethane for the preliminary extraction of these

/mixtures .....

mixtures in order to reduce the amount of soluble material which is extracted from the contaminant.

The only propellant ingredient which has been found to interfere with the procedure is lead  $\beta$ -resorcylate, which hydrolyses to produce 2:4-dihydroxybenzoic acid. This is sufficiently soluble in water to accompany resorcinol and reacts slowly with the resorcinol reagents causing somewhat high results to be obtained. Interference due to stabiliser degradation products was also considered, but this seemed unlikely in view of earlier work<sup>1,5</sup> in these laboratories which showed that the principal degradation products from resorcinol are high molecular weight compounds such as resorufin and rezazurin; and related condensation products. The formation of these products depends upon initial nitrosation of resorcinol, but there is no evidence to suggest that the intermediate nitroso compounds are themselves present in aged propellant, hence there is no reason to suspect interference from this source. However, as a precautionary measure the analysis was performed in the presence of 4-mitrosoresorcinol; no interference was observed even when relatively large proportions (ca. 30% of resorcinol content) of the nitroso compound were used.

Although the original procedure proved reasonably satisfactory for normal propellant analysis, a lack of precision in the results was often observed. During studies of the reactivity of resorcinol as a stabiliser, in which it was necessary to determine with a greater accuracy rather low rates of consumption of resorcinol after heating in solution in nitroglycerine, it became evident that variations in the concentrations and quantities of the reagents which were used could have considerable effects on the results. This might be expected from the reactions leading to structure III above which would obviously be favoured by an increase in copper ion concentration.

### 2. EXPERIMENTAL

Variations in the quantities of reagents and their effects were studied as follows:

- (a) In the original method, the methylene chloride solution obtained by extracting 1 g of propellant was washed with 3 x 20 ml of water and the resulting aqueous extract made up to 100 ml. 5 ml of this solution was mixed with 2 ml of 2% sodium nitrite solution, 20 ml of 10% acetic acid was then added, followed by 2 ml of 3% copper sulphate solution. A blank determination on the reagents was necessary to complete the estimation. It was found that the amount of sodium nitrite could be reduced to 2 ml of 0.2% solution without affecting the development of the yellow copper complex. This considerably reduced the level of the blank with a consequent improvement in the quality of the results.
- (b) The absorption spectra of solutions of the copper dinitrosoresorcinate complex as obtained by the original procedure are shown in Fig. 1. By reducing the sodium nitrite concentration and increasing the amounts of copper sulphate to a maximum of 0.5 g, the absorbance at 340 mm is increased and that at 290 mm decreased as shown in Fig. 2. The change of absorbance

/against ....

against amount of copper sulphate added is shown in Fig. 3 and it is evident that, at the low concentration required by the original method, a considerable error could result from a minor fluctuation in the amount of reagent which was added. The data showed that the method would be improved by amending the reagent concentrations as follows:

5 ml 0.2% sodium nitrite solution and 5 ml of 10% copper sulphate solution to be used for the formation of copper nitrosoresorcinate.

The detailed procedure is given in Section 3 (iii).

### 3. PROCEDURE

(i) For propellants, extract an appropriate quantity of the prepared sample for 10 hours with dichloromethane in a Soxhlet extractor. Carefully distil the extract avoiding excessive heating and redissolve in about 20 ml of chloroform. If other stabilisers, e.g. 2-NDPA, are to be determined it is convenient to make the chloroform extract up to 50 ml and pipette 25 ml for the resorcinol determination into a small separating funnel. The remaining solution is then available for further analysis.

In the preparation of samples, composite modified type propellants should be cut into small thin slices not more than 1 cm<sup>2</sup> and 1 - 2 mm thick. They should not be ground, owing to possible segregation of aluminium. Ordinary double-base propellants, and similar compositions can, of course, be ground.

For other samples or extracts dissolve an appropriate quantity 20 - 25 ml chloroform in a separating funnel.

- (ii) Extract the chloroform solution 5 times with 25 ml portions of water filtering the extracts through a cotton wool plug into a 250 ml graduated flask. Make the solution up to the mark with water.
- (iii) Pipette a suitable aliquot to contain approx. 0.5 mg of resorcinol into a 100 ml standard flask. Add 20 ml 10% acetic acid and 5 ml 0.2% sodium nitrite solution and allow to stand for 5 minutes, then add 5 ml 10% copper sulphate solution, dilute to the mark and allow to stand for 20 minutes. Measure theoptical density of the solution at 340 nm in 5 mm cells against a reagent blank prepared as above omitting the sample aliquot.
- (iv) From a calibration graph prepared with pure resorcinol calculate the percentage resorcinol in the sample as follows:

$$\%$$
 Resorcinol =  $\frac{M \times D}{W}$ 

where M = mg resorcinol in 100 ml solution

D = dilution factor
W = sample weight, g

The current calibration graph is shown in Fig. 4.

/Calibration .....

Calibration graphs produced in the absence and presence of NG by the extraction procedure given in 3 (ii) were identical.

# 4. RESULTS

The following results, expressed as weights of resorcinol in milligrammes, are from experiments in which nitroglycerine samples weighing 2.27 grammes (0.01 moles) and containing 0.5 and 2.0 molar per cent additions of resorcinol were heated at 80°C for 24 hours. The results from these experiments were self-consistent and corresponded well with the behaviour expected from such mixturesuunder these conditions.

# 0.5% Molar

Standard deviation (equal to the square root of mean variance between pairs)

$$= \sqrt{\frac{\epsilon (x_1 - x_2)^2}{2n}} = 0.04 (0.85\%)$$

# 2.0% Molar

Typical results from compatibility testing of CMCDB propellant are given below showing the stabiliser contents of propellant samples after heating along and with various Hypalon rubbers. The rates of stabiliser loss from the propellant, being no greater when heated in contact with the rubbers, show their mutual compatibility.

| Sample<br>F452/302/19     | Heating Time<br>Control<br>2 weeks 60°C<br>4 weeks 60°C | Resorcinol, % 1.18 1.20 1.21 1.08 1.10 1.08 0.96 0.96 0.96 |
|---------------------------|---|--|
| F452/302/19 + 5% Rubber 1 | 11<br>11<br>11  | 1.02<br>1.05<br>1.05<br>1.02<br>1.03                       |

### 5. CONCLUSIONS

The revised method is applicable to the determination of resorcinol in heated mixtures with nitric esters, colloidal propellants which do not contain lead  $\beta$ -resorcylate, and in compatibility studies of propellant with various

/plastics .....

plastics, rubbers and adhesives.

# 6. REFERENCES

1. Grindlay, J.W., Jeacocke, G.J.

The Stability of Colloidal Propellants: Part 6: Resorcinol as a Propellant Stabiliser

with

Grindlay, J.W., Howieson, (Miss) A.R.

Appendix: The Determination of Resorcinol ERDE Report No. 26/R/62

2. Kohlbeck, J.A.

Determination of Nitroglycerine and Resorcinol in Double-Base Propellant following Separation by Thin-Layer Chromatography. Analyt. Chem., 1965, 37, 1282.

3. Macke, G.F.

Thin Layer Chromatographic/ Spectrophotometric Analysis of Certain Components in Aged Double Base Propellant J. Chromatog., 1968, 38, 47

4. Trowell, J.M., Philpot, M.C.

Gas Chromatographic Determination of Plasticisers and Stabiliser in Composite Modified Double-Base Propellants. Analyt. Chem., 1969, 41, 166.

5. Carver, F.W.S., Gallagher, E.J.

The Application of Thin-Layer Chromatography to Propellant Analysis: Part 1: The Determination of Residual p-Nitromethylaniline in Propellants
ERDE Report No. 12/R/64

S. No. 98/70/BS






FIG.1 GRAPH OF ABSORBANCE VS WAVELENGTH
FOR COPPER DINITROSORESORCINATE
SOLUTION. ORIGINAL PROCEDURE.





FIG.2 GRAPH OF ABSORBANCE VS WAVELENGTH

FOR COPPER DINITROSORESORCINATE

SOLUTION. REVISED PROCEDURE.

CONFIDENTIAL/DISCREET

(but cleared for Canada)






FIG.4 GRAPH OF ABSORBANCE VS CONCENTRATION OF RESORCINOL.

THESE INDEX CARDS ARE SUBJECT TO THE SAME SECURITY REGULATIONS AS THE PARENT REPORT AND THEIR LOCATION SHOULD BE RECORDED ON THE INSIDE OF THE BACK COVER

| CONFIDENTIAL/DISCREET (but cleared for Canada)    | Technical Note No. 10 Explosives Research and Development Establishment THE SPECTROPHOTCHETRIC DETERMINATION OF RESORCINOL IN COLLOIDAL PROPELLANTS (C) Tunstall, F.I.H. February 1970 5 pp, no tabs, 4 figs  | photo- The method for determination of resorcinol which is based on spectrophoto- Improve metric measurement of copper nitrosoresorcinate has been modified to improve its sensitivity and reliability. | CONFIDENTIAL/DISCREET (but cleared for Canada) | CCNFICENTIAL/DISCREET (but cleared for Canada) | Technical Note No. 10 Explosives Research and Development Establishment THE SPECTROPHOTCHETRIC DETERMINATION OF RESORCINCL IN COLLOIDAL PROPELLANTS (C) Tunstall, F.I.H. February 1970 5 pp, no tabs, 4 flgs  | photo- The method for determination of resorcinol which is based on spectrophoto- metric measurement of copper nitrosoresorcinate has been modified to improve its sensitivity and reliability. | CONFIDENTIAL/DISCREET (but cleared for Canada) |
|---|---|---|--|--|---|---|--|
| CONFIDENTIAL/DISCREET<br>(but cleared for Canada) | Technical Note No. 10  Explosives Research and Development Establishment THE SPECTROPHOTCHETRIC DETERMINATION OF RESORCINOL IN COLLOIDAL PROPELLANTS (C) Tunstall, F.I.H. February 1970 5 pp, no tabs, 4 figs | The method for determination of resorcinol which is based on spectrophoto-metric measurement of copper nitrosoresorcinate has been modified to improve its sensitivity and reliability.                 | CONFIDENTIAL/DISCREET (but cleared for Canada) | CONFIDENTIAL/DISCREET (but cleared for Canada) | Technical Note No. 10  Explosives Research and Development Establishment THE SPECTROPHOTOMETRIC DETERMINATION OF RESORCINOL IN COLLOIDAL PROPELLANTS (C) Tunstall, F.I.H. February 1970 5 pp, no tabs, 4 figs | The method for determination of resorcinol which is based on spectrophotometric measurement of copper nitrosoresorcinate has been modified to improve its sensitivity and reliability.          | CONFIDENTIAL/DISCREET (but cleared for Canada) |

# DETACHABLE ABSTRACT CARDS

The Abstract Cards detached from this document are located as follows:

| 1. | Section | Signature | Date |
|----|---------|-----------|------|
| 2. | Section | Signature | Date |
| 3. | Section | Signature | Date |
| 4. | Section | Signature | Date |



CONFIDENTIAL/DISCREET (but cleared for Canada)



Information Centre Knowledge Services [dsti] Porton Down, Salisbury Wilts SP4 0.IQ Tel: 01980-613753 Fax 01980-613970

Defense Technical Information Center (DTIC) 8725 John J. Kingman Road, Suit 0944 Fort Belvoir, VA 22060-6218 U.S.A.

AD#: 399956

Date of Search: 12 December 2006

Record Summary:

Title: Resorcinol in colloidal propellants

Covering dates 1970

Availability Open Document, Open Description, Normal Closure before FOI

Act: 30 years

Former reference (Department) ERDE TN 10

Held by The National Archives, Kew

This document is now available at the National Archives, Kew, Surrey, United Kingdom.

DTIC has checked the National Archives Catalogue website (http://www.nationalarchives.gov.uk) and found the document is available and releasable to the public.

Access to UK public records is governed by statute, namely the Public Records Act, 1958, and the Public Records Act, 1967.

The document has been released under the 30 year rule.

(The vast majority of records selected for permanent preservation are made available to the public when they are 30 years old. This is commonly referred to as the 30 year rule and was established by the Public Records Act of 1967).

This document may be treated as **UNLIMITED**.